

# BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE  
TRANSPORT  
PACKAGING**



The BioTherm range of Category A and Category B solutions come complete with the Intelsius PathoPak sample transport systems and are fully compliant to IATA, ADR and DOT (49 CFR) transport regulations. The systems are qualified to ensure sample integrity for up to 96 hours against 2-8°C, 15-25°C and dry ice. Each system is provided with handles and the correct compliance labelling to enable careful handling throughout the cold chain.

## THIS PRODUCT CAN BE USED FOR

- Infectious samples
- Diagnostic samples
- Clinical trials
- Screening programmes
- Eradication schemes

## FEATURES

- Complete temperature controlled Category A and Category B compliant overpack systems
- Range of systems to transport samples at 2-8°C, 15-25°C and deep frozen to maintain sample integrity
- 4GU compliant packaging systems for use with any primary receptacle
- Ergonomically engineered, reusable, sample bottle for safe manual handling and easy extraction

For more information please contact Intelsius UK [cs@intelsius.com](mailto:cs@intelsius.com)

# BIO THERM CATEGORY A & CATEGORY B

2-8°C, 15-25°C, DRY ICE SOLUTIONS WITH PATHOPAK



**SAMPLE  
TRANSPORT  
PACKAGING**



## DRY ICE SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Empty Mass	Dry Ice Mass	Volumetric Weight L/6	Duration
BT090	BioTherm 10A Dry Ice	357 x 309 x 315 mm 14.1 x 12.2 x 12.5 in	90 dia x 170 mm 3.5 dia x 6.7 in	1 L	1.4 kg 3.1 lbs	6.9 kg 15.3 lbs	5.8 kgvol <sup>c</sup> 12.8 lbsvol <sup>c</sup>	72
BT052	BioTherm 15A Dry Ice	403 x 396 x 416 mm 15.9 x 15.6 x 16.4 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	2.4 kg 5.3 lbs	11.5 kg 25.4 lbs	11.1 kgvol <sup>c</sup> 24.4 lbsvol <sup>c</sup>	96

## HIGH PERFORMANCE 2-8°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT056P	BioTherm 30A HP 2-8 96 Summer	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	10.2 kg 22.4 lbs	17 kgvol <sup>c</sup> 37.4 lbsvol <sup>c</sup>	96
BT057P	BioTherm 30A HP 2-8 96 Winter	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	9.0 kg 19.8 lbs	17 kgvol <sup>c</sup> 37.4 lbsvol <sup>c</sup>	96

## CRT 15-25°C SOLUTIONS

Product Code	Product Name	External Dimensions	PathoPak Internal Dimensions	Payload Volume	Mass	Volumetric Weight L/6	Duration
BT105	BioTherm 30A CRT 15-25	465 x 458 x 478 mm 18.4 x 18.1 x 18.9 in	125 dia x 170 mm 4.9 dia x 6.7 in	2 L	11.9 kg 26.1 lbs	17 kgvol <sup>c</sup> 37.4 lbsvol <sup>c</sup>	96

**EUROPE:** DGP Intelsius Ltd  
1 Harrier Court, Airfield Business Park,  
Elvington, York, YO41 4EA,  
United Kingdom  
t: +44 (0) 1904 607 390  
f: +44 (0) 1904 607 391  
cs@intelsius.com

**GERMANY:** DGP Intelsius GmbH  
Eulerweg 11,  
64291 Darmstadt  
Germany  
t: +49 6151 4936650  
f: +49 (0) 6151 4936651  
info@intelsius.de

**AMERICAS:** DGP Intelsius LLC  
7696 Zionsville Road,  
Indianapolis, IN 46268,  
United States of America  
t: +1 317 452 4006  
f: +1 317 290 4621  
americasorders@intelsius.com

**IRELAND:** DGP Intelsius Ireland Ltd.  
Unit 8. Coolcaslagh Industrial Estate,  
Coolcaslagh, Killarney  
Co. Kerry, V93 K006, Ireland  
t: +353 (0) 26 45699  
salesireland@intelsius.com

**ASIA:** DGP Intelsius Sdn Bhd  
Plot 11, Lorong Bemban 2, Estet  
Perindustrian Bemban, 31000  
Batu Gajah, Malaysia  
t: +605 3666 399  
f: +605 3669 399  
compliance@intelsius.com

**INDIA:** DGP Intelsius Cold Chain Pvt. Ltd.  
Sun Srishti, C 10, Saki Vihar Road,  
Chandivali, Powai, Mumbai - 400072  
Maharashtra - India  
t: 91 22 49746536  
compliance@intelsius.com